

GABA_A Receptors

The majority of inhibitory neurotransmission in the brain is mediated by the binding of γ -aminobutyric acid (GABA) to the GABA_A receptor. These receptors mainly consist of two α subunits, two β subunits, and either a γ or δ subunit with the varying subunit compositions determining the biophysical and pharmacological properties of the receptor. Dysregulation of GABAergic neurotransmission is associated with many neurological diseases and psychiatric disorders, making these receptors attractive targets for drug discovery programs.

Developments in automated patch clamp technologies have allowed for rapid testing of large numbers of compounds against ion channel targets, however until now, the lack of an expansive panel of GABA_A cell lines has limited the identification and characterization of subtype-selective GABA_A compounds. To help overcome this, SB has developed a GABA_A receptor discovery platform comprising a comprehensive panel of 20 human GABA_A recombinant cell lines stably expressing each α subunit in combination with one of three β subunits and γ 2L. In addition this panel also includes the extra-synaptic α 4 β 3 δ and α 6 β 3 δ receptor subtypes, creating opportunities to improve subtype-specific drug efficacy, safety and tolerability.

GABA_A Receptor Discovery Platform

Utilizing the 384-well SyncroPatch (Nanon Technologies), SB has designed a comprehensive suite of flexible electrophysiology assays capable of assessing agonists, modulators and inhibitors of GABA_A receptors. This rapid, high-throughput assay format has been rigorously optimized to ensure consistent high quality data generation, supporting you throughout your discovery campaign.



Methods

HEK-293 cells stably expressing GABA_A subunit combinations were generated at SB Drug Discovery. Whole-cell patch-clamp experiments were carried out at room temperature using multi-hole chips on the SyncroPatch automated electrophysiology platform. To monitor currents, a steady state voltage pulse of -80 mV was applied in combination with the stacked addition method for rapid compound application and removal.

Using these optimized conditions and proprietary processes, GABA elicits a reproducible, concentration-dependent activation of the GABA_A receptor, with consistently comparable potencies shown across all receptor subtypes (Figure 1).




Figure 1: Example of current traces (A), time-courses (B) and EC₅₀ plots (C) for the GABA α 1 β 2 γ 2L cell line in the presence of increasing concentrations of agonist GABA (0.3-300 μ M; blue), followed by application of a saturating concentration of GABA (1 mM; grey). The normalised concentration-response curve shows a calculated value of approximately 14 μ M. Mean \pm S.E.M is shown ($n \geq 3$ cells per concentration).

Positive Allosteric Modulator Studies

In positive allosteric modulator (PAM) mode, receptor activity is initially measured in the presence of GABA (EC₁₀). GABA alone is applied three times to verify the reproducibility of the response. After pre-incubation of the compound, re-application of agonist (GABA EC₁₀) in the presence of the compound enables measurement of modulatory activity. This is followed by an application of a saturating concentration of agonist (GABA EC₁₀₀) to elicit maximal response. Using this approach, the effect of test compound can be calculated as a fold increase over GABA alone or normalised to maximum GABA response, enabling assessment of potency, efficacy and kinetics (Figure 2).

Figure 2: Effect of neurosteroid Allopregnanolone on GABA α 1 β 2 γ 2L subtype




Figure 2: Example of current traces (A) and corresponding time-courses plots (B) for the GABA α 1 β 2 γ 2L cell line in the presence of increasing concentrations of neurosteroid Allopregnanolone (0.01-1 μ M; blue), followed by application of a saturating concentration of agonist GABA (1 mM; grey). The concentration-response curve of Allopregnanolone is shown (C) as fold increase above GABA alone with a calculated EC₅₀ value of approximately 115 nM. Mean \pm S.E.M is shown ($n \geq 3$ cells per concentration). Superimposed concentration-response curves of Allopregnanolone (D) verify robust day-to-day reproducibility with EC₅₀ values within 2-fold.

SB has used its GABA_A receptor discovery platform to characterize a range of GABA_A receptor subtype-selective compounds, including α 2/3 selective anxiolytics and compounds selectively modulating α 5, α 6 and δ -subunits (Figure 3).

Figure 3: Effect of selective PAMs against GABA_A alpha subtypes




Figure 3: Normalised concentration response curves of benzodiazepine Diazepam (A), neurosteroid Ganaxolone (B), sedative-hypnotic Zolpidem (C) and null benzodiazepine agonist PZ-II-029 (D) are shown. The plots have been normalised to highlight the maximum response observed against any of the alpha subunits (α 1- α 6).

The advent of 384-well automated electrophysiology platforms has significantly increased the possibilities for GABA_A receptor drug discovery, generating a wealth of high quality electrophysiology data and guiding researchers on compound potency and subtype selection. In combination with robust, high quality cell line reagents, we have the potential to speed-up the search for lead candidates in the GABA_A receptor drug discovery process and, for the first time, offer the possibility to quickly assess the molecular mechanisms of these potential drug candidates. In the future, the aim is to achieve more selective GABA_A receptor-targeted drugs in order to improve their efficacy, safety and tolerability for a range of therapeutic purposes. As our understanding of GABA_A receptor subunit composition and molecular pharmacology improves, there is great hope for the discovery of more effective CNS therapeutics.

SB Drug Discovery's Commitment to Ion Channel Research

SB Drug Discovery is proud to play a role in these discoveries through our comprehensive suite of GABA_A receptor tools and years of experience. With state-of-the-art automated patch clamp, high-throughput screening and lead optimisation capabilities, SB Drug Discovery provides the ideal platform for in-depth analysis of GABA_A receptor modulation.

SB Drug Discovery's expertise and cutting-edge resources make us the perfect partner for collaborative research on GABA_A receptor modulation. Reach out to us today to explore the exciting possibilities of working together to advance your research goals in this promising field.



Cell Line
Generation



Drug Discovery
Screening



Ion Channel
Electrophysiology



Inflammation

SB Drug Discovery
Telford Pavilion, Todd Campus
West of Scotland Science Park
Glasgow
G20 0XA
United Kingdom

T: +44 (0)141 587 6100
enquiries@sbdrugdiscovery.com
sbdrugdiscovery.com



**SB DRUG
DISCOVERY**

A SYGNATURE DISCOVERY BUSINESS